

Installation Instructions OCT_D v120

15.09.2021

1. copy the file „DBBC3 Control OCT_D_v120.exe“ in the folder C:\DBBC\bin
2. copy the firmware-file
 - dbbc3_oct_D_1hv1_270921.bit (for DBBC3-1hv1)
 - dbbc3_oct_D_2hv2_270921.bit (for DBBC3-2hv2)into the folder C:\DBBC_CONF\FilesDBBC
3. copy the config files
 - dbbc3_config_file_oct_D_120.txt
 - oct_D_core3H_1.fila10g, oct_D_core3H_2.fila10g, ...into the folder C:\DBBC_CONF\OCT_D_120 (**new subfolder!**)
 - make sure there is a fila10g-file present for each Core3H in the system.
 - Use the example files as reference to modify your existing files from previous versions, there have been some changes in the config files from version v120.
4. copy the filter coefficient files (*_64taps.flt) into the folder C:\DBBC_CONF\OCT_D_120
5. follow the instructions in the manual
 - „Setting up the DBBC3 for OCT_D mode_v120.pdf“